# VARIXIO

## Pro Mag

EN	INSTRUCTIONS FOR USE	1
ES	INSTRUCCIONES DE USO	12
DE	GEBRAUCHSANWEISUNG	24
FR	MODE D'EMPLOI	36
IT	ISTRUZIONI PER L'USO	48
PT	INSTRUÇÕES DE UTILIZAÇÃO	59

**VARIXIO® Pro Mag**  
**INSTRUCTIONS FOR USE**

**CONTENT**

1 INTRODUCTION .....	2
2 DESCRIPTION .....	2
3 UNPACKING .....	3
4 COMMISSIONING .....	3
5 SWITCHING ON .....	4
6 INSTRUCTIONS FOR THE PREPARATION OF SCLEROSING FOAM .....	4
7 CORRECT USE .....	4
8 TECHNICAL DATA .....	5
9 EXPLANATION OF SYMBOLS AND LABELING .....	6
10 ERROR CODES.....	7
11 MAINTENANCE .....	7
12 WARRANTY AND DECLARATION OF CONFORMITY .....	9
13 INFORMATION ABOUT EMC.....	9
14 SAFETY INFORMATION .....	9
15 ELECTRICAL SAFETY REQUIREMENTS .....	10
16 ELECTROMAGNETIC SAFETEY REQUIREMENTS.....	11

## 1. INTRODUCTION

### Application and intended use

Please read this instruction manual in its entirety before using the device

### Product Description

VARIXIO® Pro Mag is an electronic, digitally programmable, magnetic stirrer, fitted with a dock in which to fit a VARIXIO® Pod capsule. The stirrer contains a rotating magnetic element that aligns with the frothing disc inside the VARIXIO® Pod capsule to make it rotate at predefined combinations of speed and time according to the gas, sclerosing agent type, and concentration parameters selected by the user.

### Models. Product references.

One model available: VARIXIO® Pro Mag, Ref. VB04ES

**Intended use:** VARIXIO® Pro Mag is intended for exclusive use with VARIXIO® Pod capsules for the preparation of foam for sclerotherapy of varicose veins.

VARIXIO® Pro Mag is an accessory device to VARIXIO® Pod capsules that must be used only by doctors or other healthcare professionals experienced in sclerotherapy protocols for the treatment of varicose veins.

**Contraindications:** No contraindications are known.


**Side effects:** No side effects are known.

## 2. DESCRIPTION

### Device Parts

VARIXIO® Pro Mag has the following components / parts:

- A Rotating mechanism: Pre-programmed magnetic rotation system with speed and time control (not visible).
- B Capsule dock: Circular compartment in which to place the VARIXIO® Pod capsule.
- C Information display: 3-line display, as an element of communication with the user.
- D Selector: Control button with rotary selection and executive function with push.
- E Power supply: External power source with 100 - 240 V AC input providing 24 V DC output.
- F Power button: Press button to start, to stop (if running) or power off the device.
- G Back button: Control button to allow navigating back to the previous selection.
- H Home button: Control button for returning to the home screen (selection Capsule Type).
- I VARIXIO® Pod capsule placed in the capsule dock.

### Cables and Accessories

- Power supply unit with interchangeable plug adapter

## 3. UNPACKING

### Unpacking


- Please unpack the device carefully
- In the case of any damage please report this immediately to your supplier or contact: [support@vbdevices.com](mailto:support@vbdevices.com).

### Delivery scope

- VARIXIO® Pro Mag stirrer unit
- Instructions for use
- Power supply unit with exchangeable plug adapters

## 4. COMMISSIONING

Pay attention to the ambient conditions listed in **8. TECHNICAL DATA**.

If these conditions are met, the unit is ready for operation once the power supply unit has been plugged in.

If the above procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged.

## 5. SWITCHING ON

Before using the appliance for the first time, carefully read sections **15. ELECTRICAL SAFETY REQUIREMENTS** and **16. ELECTROMAGNETIC SAFETY REQUIREMENTS**, which detail the essential electrical safety requirements and electromagnetic safety requirements, and **8. TECHNICAL DATA**, that provides information on location, storage and transport requirements in order to avoid malfunctioning or damage to the appliance, operator or others.

① Press “Power” button to switch on the device. The start-up screen will be displayed for several seconds. This screen shows the software version.

MAGNETIC STIRRER  
VARIXIO Pro Mag  
Version 0.0.2

② After some seconds, the capsule selection screen will appear, and the device is ready for operation. At this point, the device is ready for the user to navigate through the menu and select the appropriate parameters by turning and pressing the Selector button. **To proceed to this step see the Instructions for Use of VARIXIO® Pod.**

Capsule Type  
AIR  
O2CO2

## 6. INSTRUCTIONS FOR THE PREPARATION OF SCLEROSING FOAM

See Instructions for Use of VARIXIO® Pod.

## 7. CORRECT USE

VARIXIO® is a system for the preparation of sclerosing foam for the treatment of varicose veins. VARIXIO® consists of a customized magnetic stirring device (VARIXIO® Pro Mag) and a sterile disposable mixing capsule (VARIXIO® Pod) in which the foam is generated.

The system is designed to prepare foam with different combinations of gas types and sclerosing agents - polidocanol (POL) or sodium tetradecyl sulphate (STS) – in a range of concentrations. The gas type is determined by the gas content of the VARIXIO® Pod capsule being used. For further information please see the Instructions for Use of VARIXIO® Pod.

Note: VARIXIO® Pro Mag must be used only by doctors and healthcare professionals experienced in sclerotherapy protocols for the treatment of varicose veins.

The selection of incorrect parameters in the menu or accidental interruptions of the program for any reason may result in a sub-optimal foam, but this has not been found to represent a risk for the patient when the foam volume injected is limited to the indicated 9 ml. In any case, if for whatever reason, the foam generated shows large visible bubbles it is recommended to discard the VARIXIO Pod® capsule used and to start again with a new one.

Note: Before using the VARIXIO® Pro Mag stirring device, please also read the instructions for use of the VARIXIO® Pod capsules.

Range of use: The device is suitable only for use in interior areas of the healthcare environment.

The safety of the user cannot be guaranteed if:

- the device is operated with accessories that are not supplied or recommended by the manufacturer
- the device is operated improperly contrary to the manufacturer's specifications

## 8. TECHNICAL DATA

### Power supply unit

Input	V AC	100 ... 240
	A	1.0
	Hz	50 / 60
Output	V DC	24
		40 W.LPS (limited power source)
Protection class		2 (double insulated)



### VARIXIO® Pro Mag stirrer

Operating voltage	V DC	24
	mA	1500
Power consumption, normal operation	W	36
Power consumption, standby operation	W	2
Power output	W	30
Drive		EC flat motor
Speed range	rpm	1000-4000 (subject to program selection)
Program selection		rotating/pressing knob

Protection to EN 60529		IP 20
Contamination level		2
Overvoltage category		II
Operation at a terrestrial altitude level	m	max. 2000
Dimensions (W x D x H)	mm	122 x 178 x 54
Weight	kg	0.9
Temperature range	°C	During operation: +5 ... +40 Storage: +5 ... +40 During transport: -10 ... +60
Humidity	%	During operation: 20 ... 80 Storage: 20 ... 80 During transport: 20 ... 80

*Subject to technical changes!*

## 9. EXPLANATION SYMBOLS AND LABELING

Symbol	Explanation
	Consult instructions for use
	Follow instructions for use
	Caution
	Manufacturer
	Date of manufacturing
	Serial number
	Catalogue number
	Marking of conformity according to the European Directive
	Direct current
	WEEE (Waste from electrical and electronic equipment). Product must be collected and disposed of separately from household waste.
	Temperature limitation

### Device exterior labelling 1 - Information.

Label present on the back of the equipment next to the power connector, in accordance with the requirements established in the CE Directive

## 10. ERROR CODES

Any malfunction during operation may be identified by an error message on the display.

Proceed as follows in such cases:

- ① Switch off device using the power switch.
- ② Carry out corrective measures (see table below).
- ③ Restart the device.

Error code	Cause	Effect	Solution
<b>Error 3</b>	Inside temperature is too high	motor off	- Switch off the device and let it cool down. - Restart the device.
<b>Error 4</b>	Motor cannot run normally (cannot start or unable to reach the target speed)	motor off	- Switch off the device. - Check the structure for possible blockages.
<b>Error 47</b> <b>Error 48</b>	Overload	motor off	- Switch off the device. - Decrease the speed setting or the load.
<b>Error 50</b>	Rotor is protected to prevent overload	motor off	- Switch off the device.


If the actions described fail to resolve the fault or another error code is displayed, then take the following steps:

- Contact the service department: [support@vbdevices.com](mailto:support@vbdevices.com)
- If authorized by the service department, send the device for repair, including a short description of the fault, according to the instructions.

## 11. MAINTENANCE

### Maintenance

VARIXIO® Pro Mag is maintenance-free but will undergo normal wear.

For cleaning disconnect the main plug.

To clean the device, use only isopropyl alcohol and mild soapy water.

Use only original cables and accessories, other cables and accessories may adversely affect the EMC performance, operation and useful life of the device.

- No preventive maintenance is necessary.
- Wear protective gloves when cleaning the instruments.
- Electrical instruments must not be placed in the cleansing agent for the purpose of cleaning.
- Do not allow moisture to get into the instrument when cleaning.
- Before using a cleaning or decontamination method other than that recommended, the user must request approval from Vascular Barcelona Devices, S.L. to ensure that this method does not damage the instrument.

### Expected service life

The expected service lifetime for VARIXIO® Pro Mag is 5 years, during which it is expected to maintain its safety and essential performance.

### Repair

Please contact.

Please send in the device for repair only after it has been cleaned and is free from any contaminating materials which may constitute a health hazard.

If you require servicing, return the device in its original packaging. Please also use suitable transport packaging.

### Storage

Unplug the device for storage. Follow the steps below:

- Disconnect the cable.
- Make sure no capsule is attached.
- Store the device in its packaging in a closed and dry room.
- After prolonged storage in a cold environment, condensation may occur inside the unit when it is brought into its place of operation. Wait at least 1 hour for the device to reach room temperature and dry out before switching it on.

### Storage conditions

Temperature	+5°C – +40 °C
Humidity	Max. 80 % relative humidity
Air pressure	Min. 780 mbar

## Disposal

The device must be managed as WEEE (Waste from Electrical and Electronic Equipment), therefore it must be collected and disposed of separately from household waste.

## 12. WARRANTY AND DECLARATION OF CONFORMITY

### Warranty

In accordance with Vascular Barcelona Devices, S.L. warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. The warranty does not cover wearing parts, nor does it apply to faults resulting from improper use or insufficient care and maintenance contrary to the instructions in this operating manual.

### Declaration of conformity

Vascular Barcelona Devices, S.L. hereby declares that the VARIXIO® Pro Mag device complies with the essential requirements of the European Directives. The declaration of conformity is available upon request to: [support@vbdevices.com](mailto:support@vbdevices.com).

## 13. INFORMATION ABOUT EMC

VARIXIO® Pro Mag has been tested and found to comply with the relevant requirements for Electromagnetic Compatibility.

The use of accessories (for example, humidifiers) other than those specified for this device is not recommended. They may result in an increase of emissions or a decrease in the device's immunity.

## 14. SAFETY INFORMATION

VARIXIO® Pro Mag should not be used in any situation in which treatment with sclerosant foam is not indicated. There are no known adverse effects or complications derived from the use of VARIXIO® Pro Mag when it is used according to the operating instructions for use.

### For your protection

- Read the operating instructions in full before operation and follow the safety instructions.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- The device must be supervised at all times when in operation.

Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. There is a risk of:

- splashing liquids
- projectile parts
- body parts, hair, clothing and jewellery getting caught.

- Set up the device in a spacious area on a flat, stable, clean, non-slip, dry and fireproof surface.
- The feet of the device must be clean and undamaged.
- The device is not suitable for manual operation.
- The device may become hot when in use.
- Check the device and accessories for damage prior to each use. Do not use damaged components.

The VARIXIO® Pod capsule must always be closed when the device is in operation. Switch off the device immediately if any material leaks from the capsule and clean the device.

- The working temperature of the device must not exceed 40°C.
- Only use with VARIXIO® Pod capsules
- Do not use the device in explosive atmospheres, it is not explosion proof.
- With substances capable of forming an explosive mixture, appropriate safety measures must be applied, e.g. working under a fume hood.
- The device does not start up again automatically following a cut in the power supply.
- Safe operation is only guaranteed with the accessories described in the “Cables and Accessories” chapter.

#### For protection of the device

- The voltage stated on the name plate must correspond to the mains voltage.
- The device must only be operated with the original plug-in power supply unit.
- Protect the device and accessories from bumps and bangs.
- The device must only be opened by experts.

## 15. ELECTRICAL SAFETY REQUIREMENTS

VARIXIO® Pro Mag is designed and intended for continuous operation working with universal mains voltage of 100 – 240 V, 50 – 60 Hz.

VARIXIO® Pro Mag is protected against electric shock using grounding protection. The device is grounded through the conductor of the power cable.

**WARNING:** In order to avoid the risk of electric shock, this appliance must only be connected to a power supply with grounding protection. Power supply lines must be free from transients, current and voltage spikes, voltage drops and surges.

No panel or cover must ever be removed, except by personnel authorised by the manufacturer. VARIXIO® Pro Mag must never be left switched on, left open or unattended during maintenance of the system.

Please contact the manufacturer if for any reason the device causes tripping of the protection switch of the local power supply.

**WARNING:** If the problem persists always contact the technical service authorised by the manufacturer.

## 16. ELECTROMAGNETIC SAFETY REQUIREMENTS

VARIXIO® Pro Mag has been designed in compliance with the harmonised standard for Electromagnetic Compatibility (EMC) which states that electronic equipment must not emit electromagnetic interference that may affect radio services or the essential operation of other devices. Electronic equipment must have adequate immunity to radio frequency (RF) interference in order to offer basic safety and normal operation in the presence of electromagnetic interference.

VARIXIO® Pro Mag requires special precautions with respect to EMC and must be installed and operated in accordance with the information provided about EMC. Portable and mobile RF communication equipment may affect the VARIXIO® Pro Mag appliance during normal operation.

Only the cables and accessories supplied by the manufacturer must be used. The use of accessories, transducers and cables other than those specified may cause an increase in emissions or reduce the immunity of VARIXIO® Pro Mag.

VARIXIO® Pro Mag should not be used adjacent to or stacked on top of other equipment. If this occurs, the VARIXIO® Pro Mag appliance must be observed to check that it is operating correctly.

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Copyright © 2021 Vascular Barcelona Devices, S.L. All rights reserved.