**VARI(X)IO**

EN	INSTRUCTIONS FOR USE	4
ES	INSTRUCCIONES DE USO	9
DE	GEBRAUCHSANWEISUNG	14
FR	MODE D'EMPLOI	19
IT	ISTRUZIONI PER L'USO	24
PT	INSTRUÇÕES DE USO	29

## Description

VARIXIO® is a system for the preparation of sclerosing foam for the treatment of varicose veins. VARIXIO® consists of a customised magnetic stirring device (VARIXIO® Pro Mag) and a sterile disposable mixing capsule (VARIXIO® Pod Air) in which the foam is made.

## Presentation

VARIXIO® Pod Air are transparent plastic capsules provided sterile in single-unit bags packaged in boxes of ten units for single use. Capsules have a port, closed by a Luer stopper, for loading the sclerosing agent and removing the sclerosing foam produced. Inside the capsules there is also a magnetic rotating stirring disc. On the top of the capsule there is a removable seal for the equilibration of pressures before removal of the sclerosing foam.

VARIXIO® Pro Mag is a reusable, digitally-programmable electronic stirrer device with a docking site designed to accommodate VARIXIO® Pod Air capsules. The device contains a rotating magnetic element that aligns with the magnetic stirrer in the VARIXIO® Pod Air to make it rotate at a predefined speed and time according to the sclerosing agent type and concentration selected by the user.

Reference	Description
VB03ES	Box of 10 disposable VARIXIO® Pod Air capsules
VB04ES	Adapted magnetic stirring device VARIXIO® Pro Mag

## Indications

VARIXIO® is a system for the preparation of sclerosant foam from sclerosing agents commonly used in sclerotherapy treatment.

The concentration of sclerosant to be used for foam preparation will be decided by the doctor, taking into consideration the summary of product characteristics of the sclerosing agent and the clinical guidelines on sclerotherapy.

The system provides an alternative to the conventional method for the manual preparation of sclerosing foam, combining sterile air with different concentrations of sclerosing agents such as polidocanol (POL) or sodium tetradecyl sulphate (STS) as shown below.

Gas	Sclerosant	Concentration
Air	POL	<1%
		≥ 1%
	STS	< 0.5%
		≥ 0.5%

## Instructions for use

VARIXIO® must be used only by doctors and nurses experienced in sclerotherapy protocols for the treatment of varicose veins. The procedure must be undertaken following general principles of asepsis and handling common in sclerotherapy, following the steps indicated below:

- 1 Remove the VARIXIO® Pod Air from its bag and place it in the docking site in the VARIXIO® Pro Mag in the right or left-handed position, as preferred.

- 2 Load a 5 ml or 10 ml syringe with a standard Luer tip with 2 ml of sclerosant liquid (POL or STS). Depending on the ampoule or vial of sclerosing agent used, a hypodermic needle may be required to facilitate loading. Remove the Luer cap, attach the syringe firmly to the capsule and load the capsule with the 2 ml of sclerosant liquid.

**IMPORTANT NOTE: Silicone syringes should not be used as they negatively affect the quality and stability of the resulting foam, especially when STS is used as a sclerosing agent.**

- 3 On the VARIXIO® Pro Mag, using the dial to the right of the digital screen, select the type of capsule (always AIR when using VARIXIO® Pod Air), sclerosant (POL or STS) and concentration ( $< 1\%$  or  $\geq 1\%$  for POL, or  $< 0.5\%$  or  $\geq 0.5\%$  for STS), and press the dial to confirm the selection. Without removing the syringe from the Luer port, initiate stirring by pressing the dial again (same dial as used for the parameter selection process, to the right of the digital screen). Wait as the foam is produced.

- 4 Once stirring is completed (less than 45 seconds), the device will automatically reduce speed (without stopping completely) and show a syringe icon on the display, indicating that the foam can be withdrawn.

- 5 Before withdrawing the foam, open the top seal (A) and press in the plunger on the syringe (B) to completely empty any air from the system.

- 6 Leave the capsule on the stirrer, and as stirring continues at low speed, proceed to aspirate the foam. Withdraw up to 9 ml of foam at a speed not above 1 ml per second.

- 7 Once the desired volume of foam has been drawn, the top seal and Luer lock can be closed. The remaining foam can be kept inside the capsule for subsequent withdrawal, stirring at low speed, for up to 15 minutes.

- 8 Finally, stop the stirring by pressing on the dial and properly dispose of the capsule.

## Precautions and restrictions for use

VARIXIO® should not be used in any situation in which treatment with sclerosant foam is not indicated. There are no known adverse effects or complications derived from the use of VARIXIO® when it is used according to these instructions for use. The incorrect selection of parameters from the menu or accidental interruption to the program for any reason may result in a sub-optimal foam, but this does not represent a risk for the patient provided the foam injected is limited to the indicated 9 ml. If, for any reason, the foam generated has large visible bubbles, it is recommended to dispose of the capsule and start again.

## Expiry date and storage

VARIXIO® Pod Air are provided sterile in a bag and should not be used if this bag appears to be damaged or is not correctly sealed potentially compromising sterility. VARIXIO® Pod Air capsules are single use and should be disposed of appropriately after use. VARIXIO® Pod Air are for single use and cannot be re-sterilized once the packaging has been opened. VARIXIO® Pod Air should be stored in a dry place at temperatures not above 25°C. The manufacturer does not guarantee sterility after the expiry date.

## Symbols used

Sterile by radiation

Do not reuse

Reference Number

	Serial Number
	Lot Number
	Expiry Date
	Consult Instructions for Use
	Caution: consult documents attached
	Do not used if packaging is damaged
	Do not re-sterilise
	Storage temperature
	Manufacturer
	Manufacturing date
	European conformity
	Store in a dry place
	Recyclable product

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Copyright © 2021 Vascular Barcelona Devices, S.L.
All rights reserved