VARIXIO

Pod Air

EN	INSTRUCTIONS FOR USE	4
ES	INSTRUCCIONES DE USO	11
DE	GEBRAUCHSANWEISUNG	18
FR	MODE D'EMPLOI	25
IT	ISTRUZIONI PER L'USO	32
PT	INSTRUÇÕES DE USO	39

Description

VARIXIO® is a system for the preparation of sclerosing microfoam for the treatment of varicose veins.

VARIXIO® consists of a customised magnetic stirrer device, VARIXIO® Pro Mag, and a sterile, single use mixing capsule, VARIXIO® Pod Air, in which the microfoam is made.

Presentation

VARIXIO® Pod Air comes in sterile, single-unit bags, packed in boxes of 10 units. Sterile atmospheric gas is present inside the Pod.

The VARIXIO® Pod Air capsule is made of transparent plastic with a port, closed by a Luer cap, for the loading of the sclerosing agent and the removal of the sclerosant microfoam once formed. It contains a magnetic rotating stirrer disc and a five-micron filter.

VARIXIO® Pro Mag is a programmable and reusable digital electronic stirrer with a Pod docking site on top. The stirrer contains a rotating magnetic element that aligns with the Pod's magnetic stirrer disc to rotate it at a predefined speed and time, according to the type and concentration of sclerosing agent selected by the user.

Reference	Description
VB03ES	Box with 10 VARIXIO Pod Air
VB04ES	VARIXIO® Pro Mag programmable stirrer

Indications

VARIXIO® is a system designed to prepare sclerosant microfoam from sclerosing agents commonly used in varicose vein sclerotherapy treatments. The system represents an alternative to manual methods for generating sclerosant foam, enabling the combination of sterile atmospheric air with sclerosing agents; Polidocanol (POL) or Sodium Tetradecyl Sulphate (STS), at different concentrations, as indicated below:

Gas	Sclerosing agent	Concentration
Sterile atmospheric air	POL	≤ 1%
		> 1%
	STS	< 0.5%
		≥ 0.5%

Instructions for use

VARIXIO® must only be used by physicians or healthcare professionals with experience in performing sclerotherapy protocols for the treatment of varicose veins. The procedure must be carried out according to the general principles of asepsis and handling that are usual in sclerotherapy. The concentration of sclerosing agent to be used to prepare the microfoam will be decided by the doctor, taking into consideration the Summary of Product Characteristics of the sclerosing agent used and the clinical guidelines on sclerotherapy.

Use VARIXIO® Pod Air following the steps below:

- 1 Open the bag and remove the Pod.

Place the Pod into the dock on top of the VARIXIO® Pro Mag stirrer in the right-handed or left-handed position according to preference.

Press the power button to turn on the stirrer.

- 2 Load 2 ml of sclerosing liquid (Polidocanol [POL] or Sodium Tetradecyl Sulfate [STS]) into a Luer or Luer-Lock cone syringe, preferably with a capacity of 5 or 10 ml.

Remove the Luer cap from the access port, firmly attach the syringe to the Pod and introduce the sclerosing liquid into the Pod. You can either leave the syringe attached or detach it and cover the port with the Luer cap provided.

- 2 **NOTE:** Avoid using syringes with a high silicone content, as they negatively affect the quality and stability of the resulting microfoam, especially when STS is used as a sclerosing agent.

NOTE: DO NOT REMOVE THE TOP SEAL, it acts as a protective barrier for the Pod's interior.

- 3 Press the rotary knob to confirm the selection and turn it to the right or left to change the selection on the display. Select the **POD TYPE: Air**, the **TYPE OF SCLEROSING AGENT: POL or STS** and **CONCENTRATION**. The fourth display shows the summary of the selection. Press the rotary knob to confirm and VARIXIO Pro Mag will start the stirring process according to the selected programme.

A counter that shows the remaining microfoam preparation time will appear on the display.


To change the selection, press the arrow button to go back to previous screens or the Home button to return to the first screen.

- 4 Once the microfoam is prepared (within a maximum time of 45 seconds) the device will automatically reduce the stirring speed to a maintenance speed. A syringe icon will appear in the centre of the display, indicating that the microfoam can now be removed, and a rotating circle will appear on the left, indicating that it is still being stirred (at maintenance speed).

NOTE: To aspirate the microfoam, the device must be stirring at the maintenance speed.

- 5 Aspirate approximately 1-2 ml of microfoam with the syringe and immediately introduce it back into the capsule. Then proceed with the controlled aspiration of the microfoam to be injected, at a rate of no more than 1 ml per second. Slowly detach the microfoam-filled syringe, gently depressing the plunger to prevent air from entering the syringe.

NOTE: Remember to follow the clinical guidelines on sclerotherapy and determine the volume of microfoam to be injected based on individual risks.

- 6 Once the desired volume of microfoam has been extracted, you can stop the stirring process by pressing the rotary knob. It is recommended to aspirate the desired volume of microfoam within 15 minutes of preparation. The stirring will stop automatically after 30 minutes.

- 7 When the stirring process stops, the following message appears:

Single-use item.
Please discard the
capsule after use.

Confirm the message by pressing the rotary knob or back button and safely discard the Pod. The VARIXIO Pod Air capsule must not be used more than once.

At this point, you can choose to either switch off the VARIXIO Pro Mag or start a new stirring process with a new Pod by pressing the Start button.

NOTE: Do not reuse. Risk of infection, contamination and loss of efficacy.

- 8 VARIXIO Pro Mag will automatically switch off after 10 minutes of inactivity.

Precautions and restrictions on use

VARIXIO® Pod Air must not be used in situations for which treatment with sclerosant microfoam is not indicated. There are no known adverse effects or complications arising from the use of VARIXIO®, if used as described in these instructions for use.

VARIXIO® Pod Air is supplied sterile and must not be used if the bag containing it is damaged or not properly sealed, as its sterility may be compromised. VARIXIO® Pod Air is for single use only and cannot be re-sterilised or re-used due to risk of infection, contamination and loss of efficacy.

Selecting the incorrect parameters from the selection menu or accidentally stopping the stirrer may result in sub-optimal foam, but this does not entail any risk to the patient. However, if, for any reason, the generated foam contains large bubbles that are visible to the naked eye, it is recommended to discard the Pod and start over.

Remember to follow the clinical guidelines on sclerotherapy and to determine the volume of foam to inject based on individual risks.

VARIXIO® Pod Air must be disposed of properly after use.

Storage conditions

VARIXIO® Pod Air must be kept in a dry place at a temperature between 5 and 25°C. The manufacturer does not guarantee sterility after the expiry date.

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Copyright © 2025 Vascular Barcelona Devices, S.L.
All rights reserved